

# In Check Dial Flow Rates

Free Flow - Minimum Resistance

Diskus—30-90 l/min

Flexhaler—60-90 l/min

Autohaler—30-60 l/min

Common Low Resistance Inhaler pMDI - 30-60 l/min

## Additional Adapters

Foradil Aerolizer (F) —60-120 l/min

Handihaler (HH) —20-45 l/min

Twisthaler (THX) —60-90 l/min



Illustration1

Diskus Accuhaler is a trademark of GlaxoSmithKline  
 Flexhaler is a trademark of AstraZeneca  
 Autohaler is a trademark of 3M Healthcare  
 Aerolizer is a trademark of Novartis  
 Handihaler is a trademark of Boehringer Ingelheim  
 Twisthaler is a trademark of Merck Sharp & Dohme Corp

# In Check Dial Kit

(3579912)

Contains



(Box of 200)

D (Item Not Pictured)

Illustration2

- (A) In Check Dial (3109912)
- (B) Handihaler Resistance Adapter - HH (3109859)
- (C) Aerolizer Resistance Adapter - F (3109783)
- (D) Twisthaler Resistance Adapter- THX (not pictured) (3109860)
- (E) Inspiratory One Way Valve Mouthpiece 200/BX (3122069)
- (F) Pediatric/Small Mouth Adapter - Sample (3122220)
- (G) In Check Dial Informational CD (3570001)